

CLINICAL CHEMISTRY | HEMOGLOBINS | IMMUNODIAGNOSTICS | POINT OF CARE

Uni-Gold™ Recombigen® HIV-1/2



Rapid HIV-1/2 Test Results in 10 Minutes

BENEFITS

Uni-Gold™ Recombigen® HIV-1/2 is a highly reliable test which provides these important benefits:

- Earlier detection (IgM / IgG) for optimal disease management
- CLIA-Waived (whole blood or finger stick)
- Moderately Complex (serum or plasma)
- One year expiration from date of manufacturing



SIMPLE STEPS:

Add one drop of patient sample



Add four drops of wash solution



Read results in only 10 minutes



INTERPRETATION:

Positive Result



Negative Result



PERFORMANCE CHARACTERISTICS:

Sensitivity: HIV-1 100%; HIV-2 99.3%


Specificity: 99.7%

For full instructions, please refer to the package insert.

ORDERING INFORMATION:

Product	Part Number	Quantity	Storage
Uni-Gold™ Recombigen® HIV-1/2	1206506	20 Tests	2-27°C
Uni-Gold™ Recombigen® HIV-1/2 Controls	1206530	0.5 ml HIV 1 positive 0.5 ml HIV 2 positive 0.5 ml negative	2-8°C

CPT CODE: 86703QW



CLIA-Waived (whole blood)
Fingerstick-Approved



Distributed by:

Diagnostics Direct, LLC
One South Main Street
Cape May Court House, NJ 08210
609-368-6918
866-358-9282
diagnosticsdirect2u.com